



**Centre of Forensic Sciences  
Investigators and Submitters**

**Technical Information Sheets  
Toxicology**

November 2023

**Contents**

Contents .....	1
Introduction .....	2
Examination Strategy and Capability .....	2
Urgent Cases .....	2
Examination .....	2
Instrumentation .....	3
Interpretation .....	3
Measurement Uncertainty .....	3
Limitations .....	3
Glossary .....	4
Appendix 1 – Screening Methods .....	6
GC/MS Screen .....	6
QTOF Screen .....	8
Head-space GC-FID analysis for volatiles (screen and quantitation) .....	8
Volatile screen (qualitative only) .....	8
Appendix 2 – Drugs Requiring Targeted Analysis .....	9
Appendix 3 – Initial Analyses by Case Type <sup>a</sup> .....	10
Appendix 4 – No Method Available .....	11
Appendix 5 – Capability of Quantitative Methods .....	12
Barbiturate method (LC-MS/MS) .....	12
GHB/BHB method (LC-MS/MS) .....	12
Cannabinoid method (LC-MS/MS) .....	12
Digoxin method (LC-MS/MS) .....	12
LC-MS/MS Mix 2 .....	12
LC-MS/MS Mix 3 .....	12
LC-MS/MS Mix 4 .....	12
LC-MS/MS Mix 5 .....	13
Appendix 6 – Capability of Targeted Qualitative Methods .....	14

## Introduction

The Toxicology Section performs analyses on biological samples (e.g., blood, urine, liver) to determine the absence/presence/concentration(s) of drugs, including alcohol and poisons.

This document is intended as a convenient investigative reference but should not be relied upon as definitive or exhaustive. Please contact the Centre of Forensic Sciences (CFS) Toxicology Section for assistance with questions of an analytical or toxicological nature by e-mail or telephone (647-329-1400 or 647-329-1430). When calling please ask for the appropriate coordinator:

Coroner's Coordinator:

[CFSToxicologyCoronerCoordinator@ontario.ca](mailto:CFSToxicologyCoronerCoordinator@ontario.ca)

Criminal Coordinator:

[toxcrim@ontario.ca](mailto:toxcrim@ontario.ca)

## Examination Strategy and Capability

The screening methods employed in the Toxicology Section are:

1. Gas Chromatography/Mass Spectrometry (GC/MS)
2. Immunoassay (IA)
3. Head-Space GC analysis for volatiles
4. Quadrupole Time-of-Flight MS (QTOF)

The targeted/quantitation methods employed in the Toxicology Section are:

1. GC
2. Liquid Chromatography (LC), LC-MS/MS
3. Head-Space GC analysis for volatiles

Capabilities of screening methods are presented in Appendix 1. While these screening methods have wide-ranging capabilities, not all drugs may be reliably detected. Appendix 2 contains a list of compounds that may not be identified by the screening methods but may be detected/quantitated by targeted methods. Many of the compounds contained in this list will not be tested for unless specifically requested. If use of a specific drug is known or suspected and is relevant it should be noted in the case synopsis.

The examination strategy, i.e., determining which tests will be performed in a case, is informed by a variety of sources including case type, case history, nature of submitted samples, analytical protocols and capabilities, and discussions with clients. The initial toxicological analyses conducted for a variety of case types are presented in Appendix 3.

## Urgent Cases

Requests for expedited analyses must meet specific criteria before being accepted as an urgent case. This process requires authorization by Toxicology Section management.

## Examination

All items are visually examined on receipt to check the seal numbers (if present), the contents, and the integrity of the packaging.

## Instrumentation

### Chromatography: Gas Chromatography (GC); Liquid Chromatography (LC)

Chromatography is an analytical technique used to separate compounds based on their chemical and structural properties. GC uses a pressurized gas, while LC uses a pressurized liquid, in the separation of compounds.

**Immunoassay (IA)**

IA detects compounds in biological fluids using a reaction of an antibody or antibodies to its antigen (i.e., the drug). This technique is primarily a screening technique; however, some IA methods are semi-quantitative, e.g., acetaminophen.

**Mass Spectrometry (MS)**

MS detects, identifies, and quantitates compounds. An MS can be coupled with a GC or an LC.

**Quadrupole Time-of-Flight-MS (QTOF)**

QTOF detects and identifies compounds. A QTOF is coupled with an LC.

**Tandem MS (MS/MS)**

MS/MS detects, identifies, and quantitates compounds and is commonly coupled to a GC or LC.

**Ultraviolet and Visible (UV/VIS) Spectrophotometry**

UV/VIS spectrophotometry identifies and/or quantitates a drug based on its UV and/or visible light-absorbing properties.

**Carbon Monoxide**

Carbon monoxide is analyzed by visible spectrophotometry. Results are expressed as % carboxyhemoglobin saturation.

**Interpretation**

Quantitative results may be expressed as 1) a concentration or 2) as < or > a concentration, e.g., when sufficient for interpretation. Blood ethanol interpretations provided in reports are generally limited to cases in which the detected concentration may be associated with fatalities, may be influenced by post-mortem artefacts, may have toxic interactions with other drugs, or in the case of motor vehicle collision, associated with impairment.

**Measurement Uncertainty**

Measurements made with all scientific instruments are associated with variability. No measurement is exact but is an estimate of the true value. Calculation of measurement uncertainty (MU) employs statistical methods to determine the range of values within which the quantitative result is likely to reside. The MU provides a reasonable estimate of the variability associated with the analytical method and is based on the analysis of matrix-matched quality control samples. A minimum of 10 such analyses are used. The MU is calculated with a confidence of 95.45 per cent using a k-factor based on the degrees of freedom as determined by the Student's *t*-test and the standard deviation of the associated quality control data. The MU is expressed in the same units in which the quantitative result is reported, e.g., ng/mL, mg/L and is reported as: quantitative result  $\pm$  MU.

**Limitations**

The focus of this laboratory is drug toxicity. Clinical blood/urine chemistry analysis, e.g., electrolytes, cell counts, gas saturation, creatinine, is not performed. Analysis for antiepileptic drugs is limited to determining drug toxicity, when warranted, based on case history. This laboratory does not have validated methods to analyze some sample types, e.g., oral fluid, hair, bile, muscle, brain tissue. There are a variety of analytical issues that may prevent the detection of some of the drugs that this laboratory is commonly capable of detecting, which include:

- matrix effects
  - degree of putrefaction
  - type of sample (e.g., splenic blood)
  - post-mortem interval
  - storage conditions
- volume of sample submitted
- low concentration of the drug/sensitivity of the method

Conversely, some novel, or rarely encountered, drugs not listed in Appendix 1 may be identified by the GC/MS or QTOF screens. In this case, analytical reference material would be acquired (if available) then analysed to confirm identity. There are drugs/compounds for which the CFS Toxicology Section does not have a method, examples of which are provided in Appendix 4.

## Glossary

### Abbreviations

Analytical results are reported in terms of mg/100 mL, mg/L, or ng/mL, as shown below:

g	gram
mg	milligram
ng	nanogram
L	litre
mL	millilitre

### Breakdown Product

A compound produced either inside or outside the body that may or may not be pharmacologically active.

### Carboxyhemoglobin saturation

The percentage of hemoglobin bound by carbon monoxide.

### Central Nervous System Depression (CNS depression)

A lowering of the functional activity of the brain and/or spinal cord. Depression of the respiratory and the cardio-regulatory centres are most relevant toxicologically.

### Confirmation

The process of verifying the presence of a drug by replicate analysis using the same or different analytical technique(s). Confirmation of an immunoassay result is achieved using a more specific analytical technique.

### Coroner's Case Analytical Summary

Contains analytical results with the fatal reference and limitations. The Coroner's Case Analytical Summary is accompanied by an Interpretive Guide with information specific to this report type.

### Detected

The drug has been identified in the sample. Identification is based on criteria specific to the analytical technique.

### Fatal Reference

A minimum drug concentration at which death has been reliably reported in the forensic literature.

### Inconclusive

The presence or absence of a drug could not be determined.

### Metabolite

The product of enzymatic conversion of a drug within the body to a different compound that may or may not be pharmacologically active.

### No [other] significant findings by a [method name(s)]

This comment is inserted to provide a reference to the methods that were used. Appendices 1 and 5 can be used to identify compounds not listed and that were either not detected or the results were deemed to not be toxicologically significant, e.g., caffeine or nicotine. This may also apply to endogenous compounds, e.g., acetone < 2 mg/100 mL.

### Not Detected

The drug is either not present or is present but at an amount that cannot be discerned from other constituents in the sample.

**Post-mortem redistribution**

A phenomenon that refers to a change (either an increase or a decrease) in blood drug concentration after death; post-mortem redistribution may occur regardless of sampling site but is commonly observed as increased drug concentrations in heart blood as compared to femoral blood.

**Putrefaction**

The decomposition of organic material that involves micro-organisms.

**Report**

Contains a comprehensive summary of analytical results accompanied by interpretative conclusions.

**Tentative**

A drug has been identified by a non-specific screening method (e.g., immunoassay) but has not been confirmed. Further analysis would be required to positively identify the compound.

**Therapeutic**

The detected drug concentration is generally considered to not be toxicologically significant. The use of this term does not imply clinical efficacy.

**Traces**

The drug was detected at a concentration below that which can be reliably quantitated. The use of this term does not imply clinical efficacy.

**Unconfirmed**

A drug has been identified by a single procedure but not quantitated or confirmed by a second analysis. Unconfirmed findings may or may not be toxicologically significant.

**Appendix 1 – Screening Methods**

Drugs that can be reliably detected by screening methods

**GC/MS Screen**

<b>A</b>	
alpha-pyrrolidinovalerophenone (α-PVP)	cotinine
acetyl fentanyl <sup>2</sup>	cyclobenzaprine <sup>2</sup>
amantadine <sup>1</sup>	cyproheptadine <sup>1</sup>
amitriptyline <sup>2</sup>	<b>D</b>
amlodipine <sup>2</sup>	desipramine <sup>2</sup>
amoxapine <sup>2</sup>	dextromethorphan <sup>2</sup>
amphetamine <sup>2</sup>	dextrorphan*
amphetamine (4-fluoro)	diazepam <sup>2</sup>
anabasine	diazepam (nor) <sup>2</sup>
atomoxetine	dibucaine
atropine/hyoscyamine	dihydrocodeine
	diltiazem <sup>2</sup>
<b>B</b>	diltiazem (desacetyl) <sup>2</sup>
benzocaine	dimethyltryptamine
benzofuran (6-(2-aminopropyl), 6-APB)	diphenhydramine <sup>2</sup>
benztropine <sup>1</sup>	doxepin <sup>2</sup>
benzylpiperazine (BZP)	doxylamine <sup>2</sup>
bromo-dragonfly	<b>E</b>
brompheniramine <sup>2</sup>	ephedrine* <sup>2</sup>
bupivacaine <sup>1</sup>	estazolam
bupropion <sup>2</sup>	etizolam <sup>2</sup>
butylone/ethylone	ethylone/butylone
butyryl fentanyl	<b>F</b>
<b>C</b>	x-fluoroamphetamine
caffeine <sup>1</sup>	fluoxetine <sup>2</sup>
carbamazepine <sup>2</sup>	fluoxetine (nor) <sup>2</sup>
cathinone (cath)	flurazepam <sup>2</sup>
n-ethyl-cath	flurazepam (n-desalkyl) <sup>2</sup>
4-fluorometh-cath	fluvoxamine <sup>2</sup>
3-methoxymeth-cath	<b>H</b>
4-methyleth-cath	haloperidol <sup>1</sup>
meth-cath	hydrocodone <sup>2</sup>
chlorcyclizine	hydroxychloroquine
chlordiazepoxide <sup>2</sup>	hydroxyzine <sup>1</sup>
chloroquine	<b>I</b>
chlorpheniramine <sup>2</sup>	ibogaine
chlorpromazine <sup>1</sup>	imipramine <sup>2</sup>
cisapride	<b>K</b>
citalopram* <sup>2</sup>	ketamine <sup>2</sup>
clomipramine <sup>2</sup>	<b>L</b>
clonidine <sup>1</sup>	lamotrigine <sup>2</sup>
clozapine <sup>2</sup>	laudanose
cocaethylene	
cocaine <sup>2</sup>	
codeine <sup>2</sup>	

levamisole  
lidocaine  
loratadine  
loxapine<sup>2</sup>

**M**

maprotiline<sup>1</sup>  
meclizine<sup>1</sup>  
mefloquine<sup>1</sup>  
meperidine<sup>2</sup>  
meperidine (nor)<sup>2</sup>  
mephedrone<sup>2</sup>  
mepivacaine<sup>1</sup>  
methadone<sup>2</sup>  
methamphetamine<sup>2</sup>  
methamphetamine (4-fluoro)  
methedrone  
methotrimeprazine<sup>2</sup>  
methylenedioxyamphetamine (MDA)<sup>2</sup>  
methylenedioxyethylamphetamine (MDEA)<sup>2</sup>  
methylenedioxymethamphetamine (MDMA)<sup>2</sup>  
3,4-methylenedioxypyrovalerone (MDPV)<sup>2</sup>  
methylone<sup>2</sup>  
methylphenidate<sup>2</sup>  
metoclopramide<sup>1</sup>  
metoprolol<sup>2</sup>  
midazolam<sup>2</sup>  
mirtazapine<sup>2</sup>  
moclobemide<sup>1</sup>

**N**

nicotine<sup>1</sup>  
nortriptyline<sup>2</sup>

**O**

olanzapine<sup>2</sup>  
orphenadrine<sup>2</sup>  
oxybutynin<sup>1</sup>  
oxycodone<sup>2</sup>

**P**

paroxetine<sup>2</sup>  
pentadrone  
pentazocine<sup>2</sup>  
pentoxyphylline<sup>2</sup>  
pentylone  
phenacetin  
phencyclidine (PCP)<sup>2</sup>  
phenethylamines (2C-B, 2C-B-Fly, 2C-T-7, PEA)  
pheniramine<sup>2</sup>

phenmetrazine  
phentermine<sup>1</sup>  
piperazine, 1-3 chlorophenyl (mCPP)  
piperazine, trifluoromethylphenyl (TFMPP)  
p-fluorofentanyl  
p-methoxyamphetamine (PMA)<sup>2</sup>  
p-methoxymeth-amphetamine (PMMA)  
procaine<sup>1</sup>  
prochlorperazine<sup>2</sup>  
procyclidine<sup>1</sup>  
propoxyphene<sup>2</sup>  
propranolol<sup>2</sup>  
pseudoephedrine\*<sup>2</sup>

**Q**

quetiapine<sup>2</sup>  
quinidine<sup>1</sup>

**R**

ropivacaine

**S**

scopolamine (hyoscine)<sup>1</sup>  
sertraline<sup>2</sup>  
strychnine<sup>1</sup>

**T**

tapentadol  
terbinafine  
ticlopidine  
tramadol<sup>2</sup>  
trazodone<sup>2</sup>  
trihexphenidyl<sup>2</sup>  
trimethoprim  
trimipramine<sup>2</sup>  
triprolidine<sup>2</sup>

**V**

valeryl fentanyl  
varenicline  
venlafaxine<sup>2</sup>  
venlafaxine (O-desmethyl)<sup>2</sup>  
verapamil<sup>2</sup>

**X**

xylometazoline

**Z**

zolpidem<sup>2</sup>  
zopiclone breakdown product

**QTOF Screen**

The QTOF screen is a powerful and sensitive method that can reliably detect the drugs included in the following methods (details are listed in Appendices 5 and 6):

- LC-MS/MS Mix 2
- LC-MS/MS Mix 3 (except carfentanil)
- LC-MS/MS Mix 4
- LC-MS/MS Mix 5 (except: diflunisal, furosemide, ibuprofen, salicylate, vigabatrin)

In addition, the QTOF screen can identify psilocin. The list of drugs potentially identifiable by QTOF is too extensive to list within this document. For questions about a specific drug not listed, please contact the appropriate [case coordinator](#).

\*The GC/MS screen and QTOF screen are not capable of distinguishing racemates, therefore compounds such as dextrorphan/levorphanol, citalopram/escitalopram and ephedrine/pseudoephedrine cannot be separated. Similarly, the GC/MS screen cannot distinguish between 2-fluoroamphetamine, 3-fluoroamphetamine, and 4-fluoroamphetamine.

**Immunoassay Tests (known cross-reactivity)****Barbiturates:**

amobarbital<sup>2</sup>  
 butalbital<sup>2</sup>  
 pentobarbital<sup>2</sup>  
 phenobarbital<sup>2</sup>  
 secobarbital<sup>2</sup>

**Head-space GC-FID analysis for volatiles (screen and quantitation)**

acetone	methanol
ethanol	n-propanol (qualitative)
isopropanol	

**Volatile screen (qualitative only)**

difluoroethane	propane	acetone
dichloromethane	butane	methyl ethyl ketone
1,1,1,2-tetrafluoroethane	isobutane	isopropyl alcohol
ethyl acetate	toluene	acetaldehyde
diethyl ether	methanol	chloroform
dimethyl ether	ethanol	gasoline



**Appendix 2 – Drugs Requiring Targeted Analysis**

Compounds that may not be identified by screening methods but might be detected and/or quantitated by targeted methods.

<b>A</b>	furosemide <sup>2</sup>	<b>S</b>
antipyrine (phenazone) <sup>1</sup>	<b>I</b>	sufentanil <sup>1</sup>
<b>C</b>	ibuprofen <sup>2</sup>	<b>T</b>
carbaryl <sup>1</sup>	<b>M</b>	toluene <sup>3</sup>
carbon monoxide <sup>4</sup>	methaqualone <sup>1</sup>	<b>V</b>
cyanide <sup>2</sup>	mexiletine <sup>1</sup>	valproic acid <sup>3</sup>
<b>D</b>	<b>P</b>	vigabatrin <sup>2</sup>
diflunisal <sup>2</sup>	phenyltoloxamine <sup>1</sup>	
<b>F</b>	physostigmine <sup>1</sup>	
fenfluramine <sup>1</sup>		
formic acid <sup>3</sup>		

Methods used for the quantitation of compounds identified in the preceding appendices are denoted as follows:

<sup>1</sup> GC-NPD

<sup>2</sup> LC-MS/MS

<sup>3</sup> GC-FID

<sup>4</sup> Visible spectrophotometry

**Appendix 3 – Initial Analyses by Case Type<sup>a</sup>**

<b>Alcohol-impaired driving:</b>	Ethanol
<b>Attempted murder:</b>	dependent upon case history
<b>Confirmation of ketoacidosis:</b>	Ethanol (includes acetone), BHB
<b>Death of child &lt; 5 years of age</b>	Ethanol, QTOF Screen, LC-MS/MS Mix 3, IA cannabinoids, IA acetaminophen
<b>Drug-impaired driving:</b>	QTOF Screen, IA cannabinoids, UDM, GHB
<b>Fatal motor vehicle collision (driver) and aviation death:</b>	Ethanol, QTOF Screen, LC-MS/MS Mix 3, IA cannabinoids, CO <sup>b</sup>
<b>Fire-related death<sup>c</sup>:</b>	CO (whole blood required)
<b>Homicide:</b>	Ethanol, QTOF Screen, LC-MS/MS Mix 3, IA cannabinoids
<b>Mandatory inquest:</b>	Ethanol, QTOF Screen, LC-MS/MS Mix 3, IA cannabinoids
<b>Possible drug-related death:</b>	Ethanol, QTOF Screen, LC-MS/MS Mix 3
<b>Rule Out/exclusionary Toxicology:</b>	Ethanol, LC-MS/MS Mix 3
<b>Sexual assault<sup>a</sup>:</b>	dependent upon case history
<b>SIU death investigation:</b>	Ethanol, QTOF Screen, LC-MS/MS Mix 3, IA cannabinoids

<sup>a</sup> dependent upon sample volume

<sup>b</sup> if fire is involved

<sup>c</sup> other analyses may be performed dependent upon evidence/suspicion of intoxication

**Appendix 4 – No Method Available**

Examples of drugs/compounds for which this laboratory does not have a method

**Animal toxins**

$\alpha$ -bungarotoxin  
conotoxin  
maurotoxin  
tetrodotoxin

**Anesthetic gases**

halothane  
isoflurane  
nitrous oxide

**Curare-related toxins**

alloferine  
toxiferine  
tubocurarine

**Other**

insulin  
lead, mercury  
lithium  
polychlorinated biphenyls (PCB)  
succinylcholine  
thallium

**Appendix 5 – Capability of Quantitative Methods****Barbiturate method (LC-MS/MS)**

amobarbital (qualitative)  
 butalbital  
 pentobarbital  
 phenobarbital  
 phenytoin  
 primidone  
 secobarbital

**Cannabinoid method (LC-MS/MS)**

tetrahydrocannabinol (THC)  
 THC (11-nor-carboxy; THC-COOH)  
 THC (11-hydroxy; THC-OH, qualitative)  
 cannabidiol  
 cannabinal

**GHB/BHB method (LC-MS/MS)**

$\gamma$ -hydroxybutyrate (GHB)  
 $\beta$ -hydroxybutyrate (BHB)

**LC-MS/MS Mix 2**

benztropine  
 benzylpiperazine  
 brompheniramine  
 caffeine (semi-quantitative)  
 clonidine

ephedrine  
 haloperidol  
 ketorolac  
 loperamide (qualitative)  
 lidocaine (semi-quantitative)

mitragynine (qualitative)  
 nicotine (semi-quantitative)  
 pseudoephedrine  
 trimeprazine (qualitative)  
 warfarin

**LC-MS/MS Mix 3**

6-monoacetylmorphine (6-MAM;  
 qualitative)  
 acetyl fentanyl  
 alprazolam  
 amitriptyline  
 amphetamine  
 benzoylecgonine  
 bromazolam (qualitative)  
 bupropion  
 carfentanil  
 chlorpheniramine  
 citalopram/escitalopram  
 clonazepam  
 clonazepam (7-amino; qualitative)  
 clonazolam  
 clonazolam (8-amino; qualitative)  
 cocaethylene  
 cocaine  
 codeine  
 cyclobenzaprine  
 dextromethorphan  
 diazepam

diazepam (nor)  
 diphenhydramine  
 etizolam  
 fentanyl  
 flualprazolam  
 flubromazolam  
 flunitrazepam (7-amino)  
 fluoxetine  
 fluoxetine (nor)  
 flurazepam (n-desalkyl)  
 hydrocodone  
 hydromorphone  
 hydroxyrisperidone/paliperidone  
 isotonitazene  
 ketamine  
 ketamine (nor)  
 lorazepam  
 meperidine  
 meperidine (nor)  
 mephedrone (qualitative)  
 methadone  
 methamphetamine

methylenedioxyamphetamine  
 methylenedioxyethylamphetamine  
 methylenedioxymethamphetamine  
 midazolam  
 mirtazapine  
 morphine  
 nortriptyline  
 olanzapine  
 oxazepam  
 oxycodone  
 oxymorphone  
 paroxetine  
 pseudoephedrine  
 quetiapine  
 risperidone  
 sertraline  
 temazepam  
 tramadol (cis)  
 trazodone  
 venlafaxine  
 xylazine  
 zopiclone

**LC-MS/MS Mix 4**

alprazolam (hydroxyl)	doxylamine	naltrexone
amoxapine	duloxetine	nitrazepam
bromazepam	flunitrazepam	nitrazepam (7-amino)
buprenorphine	flunitrazepam (N-desmethyl)	orphenadrine
butyryl fentanyl	flurazepam	PCP
chlordiazepoxide	fluvoxamine	pentazocine
chlorpromazine	furanyl fentanyl	pheniramine
clobazam	imipramine	promethazine
clomipramine	levorphanol/dextrorphan	propoxyphene
clozapine	(qualitative)	triazolam
demoxepam	loxapine	triazolam (hydroxy)
desipramine	methylenedioxypropylvalerone	trimipramine
desomorphine	methotrimeprazine	U-47700
diltiazem	methylone	venlafaxine (O-desmethyl)
diltiazem (desacetyl)	methylphenidate	ziprasidone
doxepin	naloxone	zolpidem

**LC-MS/MS Mix 5**

acebutolol	gabapentin	prochlorperazine
acetaminophen	guaifenesin	propafenone
amiodarone	ibuprofen	propranolol
amlodipine	labetalol	pseudoephedrine
atenolol	lamotrigine	salicylate
baclofen	methocarbamol	topiramate
carbamazepine (qualitative)	metoprolol	verapamil
diflunisal	naproxen	vigabatrin
furosemide	pregabalin	

**Appendix 6 – Capability of Targeted Qualitative Methods****Urine Drug Mix (UDM; LC-MS/MS)**

6-monoacetylmorphine (6-MAM)	etizolam	morphine
acetyl fentanyl	fentanyl	morphine-3-glucuronide
acetyl norfentanyl	fentanyl (nor)	morphine-6-glucuronide
alprazolam	flualprazolam	naloxone
amitriptyline	flubromazolam	naltrexone
amlodipine	flunitrazepam	nitrazepam
amoxapine	flunitrazepam (7-amino)	nitrazepam (7-amino)
amphetamine	flunitrazepam (N-desmethyl)	nortriptyline
baclofen	fluoxetine	olanzapine
benzoylecgonine	fluoxetine (nor)	orphenadrine
bromazepam	flurazepam	oxazepam
brompheniramine	flurazepam (n-desalkyl)	oxazepam glucuronide
buprenorphine	fluvoxamine	oxycodone
buprenorphine glucuronide	furanyl fentanyl	oxymorphone
bupropion	gabapentin	paroxetine
butyryl fentanyl	GHB	pentazocine
carfentanil	heroin	phenazepam
chlordiazepoxide	hydrocodone	phencyclidine
chlorpheniramine	hydromorphone	pheniramine
citalopram/escitalopram	hydromorphone-3-glucuronide	pregabalin
clobazam	hydroxyalprazolam	propoxyphene
clomipramine	hydroxytriazolam	propranolol
clonazepam	imipramine	pseudoephedrine
clonazepam (7-amino)	ketamine	quetiapine
clonazolam	ketamine (nor)	risperidone
clonazolam (8-amino)	lamotrigine	sertraline
clozapine	levorphanol/dextrorphan	tapentadol
cocaethylene	lidocaine	temazepam
cocaine	lorazepam	temazepam glucuronide
codeine	lorazepam glucuronide	THC-COOH
codeine-6-glucuronide	loxapine	THC-COOH glucuronide
cyclobenzaprine	mepidine	topiramate
demoxepam	mepidine (nor)	tramadol (cis)
desipramine	mephedrone	trazodone
desomorphine	methadone	triazolam
dextromethorphan	methamphetamine	trimipramine
diazepam	methylenedioxyamphetamine	U-47700
diazepam (nor)	methylenedioxyethylamphetamine	venlafaxine
diltiazem	methylenedioxymethamphetamine	venlafaxine (O-desmethyl)
diltiazem (desacetyl)	methylenedioxypropylvalerone	zaleplon
diphenhydramine	methylone	ziprasidone
doxepin	methylphenidate	zolpidem
doxylamine	metoprolol	zopiclone
duloxetine	midazolam	
ephedrine	mirtazapine	